



Sterisart® NF – 16478

Adapter for Medical Devices with Luer Lock Connector



International pharmacopeias require the complete sterility of pharmaceutical products that are injected into the blood stream or that otherwise enter the body below the skin surface. Manufacturers of such products are required to supply proof of the sterility of the final product batch. The preferred method for sterility testing is the membrane filter method – any microorganism present is concentrated on the membrane filter surfaces in the sample containers. In the final step nutrient media are added to these containers and the containers are incubated under the conditions specified.

The Sterisart® NF system 16478 is equipped with a dual-needle metal spike with a vent filter that is connected to the Sterisart® NF system via a single hose and a Y-distributor. A second shorter hose is attached to the other end of the Y-distributor fitted with a female luer lock connector. Using this connector, it is possible to draw samples from collapsible medical devices like bags fitted with luer lock male connectors. The dual-needle metal spike with vent filter is used to transfer rinsing solutions and/or nutrient media from closed vials.

The Sterisart® NF type 16478 is available in a special isolator version with the model suffix "gamma", which is indicated by the last letters of the order number: GBD.

Ideally the systems are used in combination with the Sterisart® Universal pumps 16419 oder 16420. These can be used in clean rooms and clean benches or can be installed in isolators.

Sterisart® NF systems can also be used with other manufacturers' pumps which have been specially designed for disposable closed sterility testing systems. This requires an adapters to ensure the correct position of Sartorius NF units in the pumps of other manufacturers: 16412V (for older pumps) or 1ZG---0014 (for the latest pump versions).

Specifications

Pore size of the Sartochem® membrane filter	0.45 µm, tested with <i>Serratia marcescens</i>
Filter area	15.7 cm ² in each Sterisart® container
Flow rate (for water)	500 ml/min at 1 bar (approx. 15 psi)
Pore size of the air filters	0.2 µm PTFE, validated acc. to HIMA for the retention of <i>B. diminuta</i>
Sample container capacity	120 ml (graduation marks at 50, 75 and 100 ml)
Max. operating pressure	3 bar (approx. 44 psi) at 20 °C
Max. operating temperature	50 °C
Sterilization	Gamma irradiation

По вопросам продаж и поддержки обращайтесь:

Астана+7(7172)727-132, Волгоград(844)278-03-48, Воронеж(473)204-51-73, Екатеринбург(343)384-55-89,
Казань(843)206-01-48, Краснодар(861)203-40-90, Красноярск(391)204-63-61, Москва(495)268-04-70,
Нижний Новгород(831)429-08-12, Новосибирск(383)227-86-73, Ростов-на-Дону(863)308-18-15, Самара(846)206-03-16,
Санкт-Петербург(812)309-46-40, Саратов(845)249-38-78, Уфа(347)229-48-12

sst@nt-rt.ru || sartorius.nt-rt.ru

Chemical compatibility of the components (24 hour contact at 20 °C)

	Sartochem® membrane	Tubing	Sterisart® container
Aqueous solutions	compatible	compatible	compatible
Alcohols	compatible	compatible	compatible
Aliphatic hydrocarbons	compatible	compatible	compatible
Oils	compatible	compatible	compatible
Salt solutions	compatible	compatible	compatible
Ether	compatible	compatible	compatible
Weak acids	compatible	compatible	limited compatibility
Weak bases	limited compatibility	compatible	limited compatibility
Aromatic hydrocarbons	compatible	not compatible	not compatible
Esters	compatible	not compatible	limited compatibility
Halogenated hydrocarbons	compatible	not compatible	not compatible

Ordering information

Name	Description	Quantity	Order number
Sterisart® NF gamma	Adapter for medical devices with luer lock male connectors (individually gamma sterilized, double-packed, optimal for use in isolators)	Box of 10	16478 GBD
Sterisart® Universal pump	Basic version: Peristaltic pump for use with sample containers in sterility testing	1	16419
Sterisart® Universal pump	Upgraded version with display user software	1	16420
Pump adapters	For use with Sterisart® systems in older competitor pumps. Please clarify pump type used.	2	16412 V
Pump adapters	For use with Sterisart® systems in newer competitor pumps. Please clarify pump type used.	1	1ZG---0014

Additional Sterisart® units and accessories are available upon request.

Comprehensive validation support literature is available on request.
Publ. No. SL-4507-e

По вопросам продаж и поддержки обращайтесь:

Астана+7(7172)727-132, Волгоград(844)278-03-48, Воронеж(473)204-51-73, Екатеринбург(343)384-55-89,
Казань(843)206-01-48, Краснодар(861)203-40-90, Красноярск(391)204-63-61, Москва(495)268-04-70,
Нижний Новгород(831)429-08-12, Новосибирск(383)227-86-73, Ростов-на-Дону(863)308-18-15, Самара(846)206-03-16,
Санкт-Петербург(812)309-46-40, Саратов(845)249-38-78, Уфа(347)229-48-12

sst@nt-rt.ru || sartorius.nt-rt.ru