



sartorius stedim  
biotech

## MidiCaps<sup>®</sup> and Capsules – 0.2 μm

The New Sterilizing Grade Filters

Membrane Filter Cartridges



# Sartopore<sup>®</sup>

platinum

A New Vision of Perfection.  
A New Class of Sterile Filtration.

**По вопросам продаж и поддержки обращайтесь:**

Астана+7(7172)727-132, Волгоград(844)278-03-48, Воронеж(473)204-51-73, Екатеринбург(343)384-55-89,  
Казань(843)206-01-48, Краснодар(861)203-40-90, Красноярск(391)204-63-61, Москва(495)268-04-70,  
Нижний Новгород(831)429-08-12, Новосибирск(383)227-86-73, Ростов-на-Дону(863)308-18-15, Самара(846)206-03-16,  
Санкт-Петербург(812)309-46-40, Саратов(845)249-38-78, Уфа(347)229-48-12

**sst@nt-rt.ru || sartorius.nt-rt.ru**

# Description

Sartopore® Platinum sterilizing grade MidiCaps® and Capsules are self-contained ready-to-use, sterile filter units for sterilizing-grade filtration in the pharmaceutical | biotech industry. These filters contain a unique heterogeneous double layer of hydrophilized polyethersulfone membranes and are specially designed to meet the requirements for filtration of a broad range of pharmaceutical products.

The complete new and innovative technologies which are incorporated in these filters lead to outstanding and unique performance data. Using Sartopore® Platinum the critical step of sterilizing-grade filtration will reach a yet unmatched quality, performance and cost efficiency:

- outstanding total throughput and permeability
- unmet wettability
- broad chemical compatibility (pH 1-14)
- high thermal resistance
- exhibiting very low protein binding
- low extractables level
- reliable integrity testing

## New Surface Modification

A new and patented membrane hydrophilization process is used to permanently modify the membrane surface. In this process, a thermally exceptionally stable and hydrophilic polymer is directly grafted to the inner and outer surface of the membrane.

This new technology provides those membrane surface properties that are responsible for the outstanding wettability and low protein binding character of the Sartopore® Platinum membrane, even after thermal and chemical stress.

## Excellent Wettability

Sartopore® Platinum MidiCaps® and Capsules can be easily wetted with lowest flushing volume needed. For reliable wetting you only need less than 1 liter for a size 9 MidiCaps® (0.26 m<sup>2</sup>).

This excellent wetting behavior will thus directly reduce the filtration costs of your process.

## Applications

Typical applications include sterilizing grade filtration of:

- Therapeutics
- Biological Fluids
- Injectables
- Media
- Buffers
- Chemicals
- Cleaning and sanitizing agents

## Reliable Integrity Testing

Imperfect wetting is the most frequent reason for failed integrity tests. In such cases filters have to be reliably re-wetted and tested again. Besides the cost factor for the additional test also a certain risk is present to lose a complete product batch when even repeated tests fail again.

The extraordinary wetting behavior of Sartopore® Platinum thus avoids this important risk factor. Using Sartopore® Platinum the process of integrity testing meets a yet unmet degree of reliability.

## Compatibility

The polyethersulfone membrane can be used for the filtration of liquids with a pH range between 1 and 14. It is not affected even by multiple sterilizing cycles. Therefore, Sartopore® Platinum MidiCaps® and Capsules are suitable for most filtration applications in pharmaceutical and biotechnological processes.

## Scalability

The flow rates and total throughput per squaremeter membrane area are the same for the different filter sizes. Therefore, Sartopore® Platinum can be perfectly up-scaled from R&D to large process scale.

## Microbial Retention

Sartopore® Platinum MidiCaps® and Capsules are fully validated as sterilizing grade filter elements according to HIMA and ASTM F-838-05 guidelines.

## Packaging

All MidiCaps® and Capsules are sterile and double wrapped for easy and safe clean-room transfer.

## Quality Control

Each individual filter element is integrity tested by diffusion and bubble point prior to release, assuring highest quality and absolute reliability.

## Documentation

Sartopore® Platinum MidiCaps® and Capsules are designed, developed and manufactured in accordance with an ISO 9001 certified Quality Management System. A Validation Guide and Extractables Guide are available for compliance with regulatory requirements.



# Specifications

## Materials of Construction

Prefilter membrane	Polyethersulfone, asymmetric
Endfilter membrane	Polyethersulfone, asymmetric
Support fleece	Polypropylene
Core	Polypropylene
End Caps	Polypropylene
Capsule Housing	Polypropylene
O-rings	Silicone

## Pore Size

0.45 µm + 0.2 µm (double layer)

## Available Sizes | Filtration Area

### Capsules

Size 4 0.021 m<sup>2</sup> | 0.226 ft<sup>2</sup>

### MidiCaps®

Size 7 0.065 m<sup>2</sup> | 0.67 ft<sup>2</sup>

Size 8 0.13 m<sup>2</sup> | 1.4 ft<sup>2</sup>

Size 9 0.26 m<sup>2</sup> | 2.8 ft<sup>2</sup>

Size 0 0.52 m<sup>2</sup> | 5.6 ft<sup>2</sup>

## Available Connectors Capsule Size 4

SS, SO, OO, NO

S: TC Flange 25 mm (3/4")

O: multiple stepped hose barb

N: 1/4" NPT threat

## Available Connectors MidiCaps®

SS, SO, OO, FF, FO, HH

S: TC Flange 50 mm (1 1/2")

O: 1/2" Single-stepped hose barb

F: TC Flange 25 mm (3/4")

H: 1/4" Multiple stepped hose barb

## Operating Parameters Capsule Size 4

Max. allowable differential pressure 4 bar | 58 psi at 20°C  
2 bar | 29 psi at 50°C

Max. allowable back pressure 2 bar | 29 psi at 20°C  
1.5 bar | 21.8 psi at 50°C

## Operating Parameters MidiCaps®

Max. allowable differential pressure 5 bar | 72.5 psi at 20°C  
3 bar | 29 psi at 50°C

Max. allowable back pressure 2 bar | 29 psi at 20°C  
1.5 bar | 21.8 psi at 50°C

## Extractables

Sartopore® Platinum 0.2 µm rated MidiCaps® and Capsules meet, or exceed the requirements for WFI quality standards set by the current USP. A detailed Extractables Guide is available on request.

## Regulatory Compliance

100% individually integrity tested (Diffusion Test and Bubble Point Test).

Integrity test correlated with HIMA/ASTM F838-05 Bacterial Challenge Test.

Non pyrogenic according to USP Bacterial Endotoxins.

Passes USP Plastics Class VI Test.

Non-fibre releasing according to 21 CFR.

## Sterilization

### Autoclaving

134°C, 2 bar | 29 psi, 30 min

No in-line steam sterilization

### Sterilization Cycles


Autoclaving: 25 (MidiCaps®)  
3 (Capsules)

## Technical References


### Validation Guide

SPK5795-e

Water Flow Rates for MidiCaps®



Water Flow Rates for Capsules Size 4



# Ordering Information

Order Code	Pore Size [ $\mu\text{m}$ ]	Pack Size [pieces]	Test Pressure [bar   psi]	Max. Diffusion [ml/min]	Min. Bubble Point [bar   psi]
<b>Capsules Size 4</b>					
5491307H4--**--B	0.2 $\mu\text{m}$	5	2.5   36.25	1.1	3.5   50.7
<b>MidiCaps®</b>					
5495307H7--**--A	0.2 $\mu\text{m}$	4	2.5   36.25	4	3.5   50.7
5495307H8--**--A	0.2 $\mu\text{m}$	4	2.5   36.25	5	3.5   50.7
5495307H9--**--A	0.2 $\mu\text{m}$	4	2.5   36.25	7	3.5   50.7
5495307H0--**--V	0.2 $\mu\text{m}$	2	2.5   36.25	14	3.5   50.7

\*\* Connector type

**По вопросам продаж и поддержки обращайтесь:**

Астана+7(7172)727-132, Волгоград(844)278-03-48, Воронеж(473)204-51-73, Екатеринбург(343)384-55-89,  
Казань(843)206-01-48, Краснодар(861)203-40-90, Красноярск(391)204-63-61, Москва(495)268-04-70,  
Нижний Новгород(831)429-08-12, Новосибирск(383)227-86-73, Ростов-на-Дону(863)308-18-15, Самара(846)206-03-16,  
Санкт-Петербург(812)309-46-40, Саратов(845)249-38-78, Уфа(347)229-48-12

**sst@nt-rt.ru || sartorius.nt-rt.ru**