



**sartorius stedim**  
biotech

## Sartocheck® mini

Filter Integrity Tester for Food & Beverage Applications



### Description

The automatic filter integrity tester Sartocheck® mini can be used to verify the integrity of membrane filters which are used in the food & beverage environment.

Taken those specific needs into account, this unit offers the following main features:

- Automatic filter integrity tester
- Pressure Drop Test &
- Diffusion Test
- Small, portable unit
- 19 different test programs
- 100 test results to be stored
- LCD display
- Automatic venting after the test
- Thermo-printer (57 mm paper)
- Easy and reliable data transfer to PC
- High capacity batteries for up to 4 hrs work
- For filter housings with a net volume between 0.1 and 999 L
- Protection rating IP50
- Incl. bag and case

### Test Result Documentation

Test results are automatically printed using the built-in thermo printer. An additional port allows the connection of an external printer.

### Data Storage

The unit stores up to 100 test results in the internal memory. To avoid the oldest data to be overwritten, electronic data can be stored on a connected PC with user-friendly software. The same software can be used for programming the device.

### Technical Specifications

Power requirements	100–240 V AC, 50 60 Hz
Max. Power Input	20 W
Max. inlet pressure	4500 mbar
Dimensions	315×150×280 mm
Weight	ca. 3.900 g
Languages	German, English, French, Italian, Spanish, Portuguese

### Operating Conditions

Temperature	3–30°C.
Humidity	5–95% rel.

### Measuring Ranges

Test pressure	0–3900 mbar
Max. inlet pressure	0–4500 mbar
Net volume	0.1–999 L

### Measuring Accuracy

Rel. deviation pressure measurement	< 0.2%
Abs. deviation pressure measurement	max. ±4 mbar (@20°C)

### Interfaces

External printer	Centronics 25 pol
Communication port	232, 9 pole male

### Equipment Supplied

- Sartocheck® mini integrity test unit
- Low volume adapter for net volumes <5 L
- Printer paper (4 rolls)
- Pressure inlet tubing (18104)
- Pressure outlet tubing (18103)
- Carrier bag (soft case)
- Hard case

### Order Information

26292---01

### По вопросам продаж и поддержки обращайтесь:

Астана+7(7172)727-132, Волгоград(844)278-03-48, Воронеж(473)204-51-73, Екатеринбург(343)384-55-89,  
Казань(843)206-01-48, Краснодар(861)203-40-90, Красноярск(391)204-63-61, Москва(495)268-04-70,  
Нижний Новгород(831)429-08-12, Новосибирск(383)227-86-73, Ростов-на-Дону(863)308-18-15, Самара(846)206-03-16,  
Санкт-Петербург(812)309-46-40, Саратов(845)249-38-78, Уфа(347)229-48-12

[sst@nt-rt.ru](mailto:sst@nt-rt.ru) || [sartorius.nt-rt.ru](http://sartorius.nt-rt.ru)